

Product Description

PICC-Nate® catheters and introducers are available in a variety of kit configurations, in “Basic” or “Complete” pre-assembled procedure trays.

Basic Catheter Kit Contents

1 Single Lumen PICC Catheter 1 Peel-Away Introducer

Complete Catheter Tray Contents:

1 Single Lumen PICC Catheter	2 each: sharps disposal devices
1 Peel-Away Introducer	1 each: Scissors
2 each: Tape measure	1 each: Syringe with 18 Ga Needle
1 pkg: Povidone Iodine Swabsticks (3 per pack)	3 each: Tape Strips
8 each: 2”x2” Gauze Pads	1 each: Instructions For Use
8 each: 4”x4” Gauze Pads	1 each: Poly-Lined Drapes
1 Fenestrated Drape with Tape	1 Transparent Dressing
1 each: Tourniquet	1 each: Micro Forceps (Half Curved Fine Point)

Ordering Information and Specifications:


1.1 Fr (28 Ga/23cm)	Basic Kit (10/bx)	Complete Tray (10/bx)	Cath Only (10/bx)	Approx. Prime (mL)	
TecoFlex®	PB-1PP	-	-	0.11	
1.9 Fr (24 Ga/35cm)	Basic Kit (10/bx)	Complete Tray (10/bx)	Cath Only (10/bx)	Approx. Prime (mL)	
Silicone	PB-2PS	PC-2PS	P-2S	0.07	
TecoFlex®	PB-2PP	PC-2PP	P-2P	0.12	
Silicone Ext Hub	PB-2LPS	PC-2LPS	P-2LS	0.11	
3.0 Fr (20 Ga/35cm)	Basic Kit (10/bx)	Complete Tray (10/bx)	Cath Only (10/bx)	Approx. Prime (mL)	
Silicone	PB-3PS	PC-3PS	-	0.12	
TecoFlex®	PB-3PP	PC-3PP	-	0.16	
Silicone Ext Hub	PB-3LPS	PC-3LPS	-	0.17	
Introducers			Catalog Number	Gauge	Qty
1.1 Fr			PEEL1	26	5 / box
1.9 Fr			PEEL2	22	5 / box
3.0 Fr			PEEL3	19	5 / box
Catheter Repair Kits			Catalog Number		Qty
1.9 Fr			P-R2		10 / box
3.0 Fr			P-R3		10 / box



Utah Medical Products Inc.
 7043 South 300 West
 Midvale, Utah 84047
 (801) 566-1200
 Toll Free 1-800-533-4984
 customer@utahmed.com
 www.utahmed.com

PICC-Nate[®]

Peripherally Inserted Central Catheter Family
1.1 Fr, 1.9 Fr, 3.0 Fr



Specifically designed to minimize trauma to the neonatal patient.

- * Slender, low-profile hub reduces patient discomfort.
- * Soft, strong silicone shaft material provides long-term friendly care.
- * Suture wing is profiled for excellent hub securement.
- * Catheter flow rates provide excellent infusion therapy.
- * Accurate catheter positioning accomplished with depth markings at every centimeter.
- * Reduced catheter shaft diameters minimize puncture wound.
- * Barium sulfate-loaded catheter body provides radiopacity for placement visualization.



- * Double-beveled needle tips allow precise vessel entry.
- * Funneled entry guides PICC introduction.

Utah Medical Products, Inc. provides a broad range of critical care devices including the Gesco[™] family of neonatal products.



UTAH MEDICAL
PRODUCTS INC.



All Gesco[™] products are Latex-free