

# Dialy-Nate®

Complete Family of Peritoneal Dialysis Sets  
For Infants



Safer and more effective than piecing together separate tubing and components

Closed system greatly reduces risk of peritonitis

Saves nursing time and worry about researching hospital protocols

- \* Preassembled, closed sterile system
- \* Graduated (150ml) administration burette with injection site and air vent
- \* Three bag spikes for mixing infusate
- \* Bonded helical warming coil for convenient warming w/heating device (e.g., K-Pad, heat lamp, etc.)
- \* In-line three-way flow control stopcock
- \* Flashbulb for drainage sampling
- \* Finely graduated dialysate meter with overflow bag

## ORDERING INFORMATION

Dialy-Nate with bonded helical warming coil  
Product No: 4000507 (priming vol: ~63 ml)  
Packaging: 1 per box



UTAH MEDICAL  
PRODUCTS INC.




UNITED STATES  
7043 South 300 West  
Midvale, Utah 84047  
(801) 566-1200  
Toll Free 1-800-533-4984  
Fax (801) 566-2062  
www.utahmed.com

EUROPE  
Utah Medical Products, Ltd.  
Athlone Business & Technology Park  
Athlone, County Westmeath  
Ireland  
Tel: +353-90-647-3932

# Dialy-Nate®

With Expandable Warming Coil

The Only Pre-assembled, Closed Peritoneal Dialysis System  
for Neonatal/Pediatric Patients



## ORDERING INFORMATION

Dialy-Nate with unbound warming coil

Product No: 4000517 (priming vol: ~48 ml)

Packaging: 1 per box



UTAH MEDICAL  
PRODUCTS INC.



UNITED STATES  
7043 South 300 West  
Midvale, Utah 84047  
(801) 566-1200  
Toll Free 1-800-533-4984  
Fax (801) 566-2062  
www.utahmed.com

EUROPE  
Utah Medical Products, Ltd.  
Athlone Business & Technology Park  
Athlone, County Westmeath  
Ireland  
Tel: +353-90-647-3932