


One Gentle Insertion!

30 DAYS...




PRESORT STANDARD
U. S. POSTAGE
PAID
PERMIT NO. 3280
SALT LAKE CITY, UT

**UTAH MEDICAL PRODUCTS INC.**
7043 South 300 West ♦ Midvale, Utah 84047
800-533-4984 ♦ 801-566-1200
www.utahmed.com



Nutri-Cath 5 Fr, 6.5 Fr, and 8 Fr catheter sizes have both an end port at the rounded tip, and offset, opposing side ports to provide accurate positional delivery of enteral solutions and obstruction-free sampling.

Tip Detail



NUTRI-CATH®:

A DEVELOPMENTALLY FRIENDLY silicone long-term enteral feeding catheter from Gesco™ and Utah Medical Products, Inc.

Because providing for a neonate's most basic need shouldn't be traumatic.

Nutri-Cath can indwell for up to 30 days without adverse effects. It avoids multiple placements with a non-stiffening, minimally clogging, biocompatible silicone material. Multiple placements are highly traumatic to a baby's delicate tissues and can lead to oral aversion and hyperactive gag reflexes.

- Non-wettable silicone material will not perforate or adhere to lining of stomach during sampling.
- Less interruption of important sleep cycles because it is softer.
- Barium sulfate-loaded catheter is highly visible under x-ray for easy verification of tip placement.
- Depth markings graduated every centimeter.
- Orange hub identifies enteral-only connection to prevent accidental medication infusion.

Choice of standard luer lock connector for continuous feeding or oral dose connector for intermittent feeding.



Nutri-Cath Silicone Enteral Feeding Catheters

Packaged 10/bx

Standard Hub	3.5 Fr	5.0 Fr	6.5 Fr	8.0 Fr
40 cm	4163519	4165019	4166519	4168019
58 cm	4163529	4165029	4166529	4168029
93 cm	4163539	4165039	4166539	4168039

Nutri-Lok® Hub	3.5 Fr	5.0 Fr	6.5 Fr	8.0 Fr
40 cm	4153517	4155017	4156517	4158017
58 cm	4153527	4155027	4156527	4158027
93 cm	4153537	4155037	4156537	4158037

For more information or to place an order, visit www.utahmed.com/nutricath.htm

or call:
1-800-533-4984

DEHP FREE