

Accessories for Improved Outcomes



Arom-Cot®, Patent 6,027,511

Arom-Cot®: Non-latex digital amniotomy device that does not compromise patient confidence and relaxation.

- Optimal control: the finger is flexible allowing more precise hook placement
- Effectively ruptures membranes in any patient position
- Emulates a vaginal exam



Muc-X™ Aspirator

MUC-X™ Aspirator: Newborn Oropharyngeal and nasal aspirator that easily, safely and reliably aspirates meconium or bodily fluids without clinician exposure.

- Provides a smooth, beveled tip that is gentle to infant
- Allows for one-handed operation

Abcorp™ FM Belts & Straps: Latex-free fetal monitoring abdominal belts and stockinettes that offer a secure yet comfortable fit.

- Broad selection of buttonhole, plush elastic, soft strap, foam and mesh designs
- Stockinettes available in precut or in a roll
- Also provide pink, blue baby beanies and foam circumcision straps



Abcorp™ Belts & Straps



Quick Connect Plus™

Quick Connect Plus™: FSE for fetal heart monitoring that meets FDA safety standards for lead wires.

- Single helix configuration yields secure attachment
- Locking mechanism that secures the electrode during attachment
- Specially designed leg plate attachment pads secure electrode cables and minimize patient discomfort



CordGuard®, Patents 5,190,556 and 5,575,796

CordGuard®: A unified system designed for umbilical cord clamping, cutting, and blood collection.

- Designed to reduce the risk of needle sticks and clinician exposure to blood
- Assists in placental delivery
- Ensures a clean neonatal blood sample

UTAH MEDICAL
PRODUCTS INC.



For more information or
to place an order, visit:

www.utahmed.com

or call:

1-800-533-4984

- Safe delivery...



UTAH MEDICAL
PRODUCTS INC.

7043 South 300 West ♦ Midvale, Utah 84047
800-533-4984 ♦ 801-566-1200
www.utahmed.com

PRESORT STANDARD
U. S. POSTAGE
PAID
PERMIT NO. 3280
SALT LAKE CITY, UT

ph 56268 rev 082707

Effective products

