

# Uri-Cath™

## Closed, Pre-assembled Urinary Drainage System

Because babies in the NICU *require*  
gentle urinary drainage

Accurate measurement of urinary output  
for neonates and pediatric patients

Soft silicone material minimizes tissue  
irritation and encrustation



- Pre-assembled, sterile system reduces risk of infection and reduces valuable nursing time
- Stopcock with sampling access site allows for non-contaminated midstream specimen sampling
- 150 ml burette for accurate urine output measurement



Utah Medical Products, Inc manufactures a broad range of neonatal critical care devices including the Gesco™ family of products

UTAH MEDICAL  
PRODUCTS INC.





## Uri-Cath™

### Product Description

The Uri-Cath™ tray consists of a pre-connected silicone urinary drainage catheter, three-way stopcock for sampling purposes, and a 150 ml burette for accurate urine output measurement.

### Design Features

- Biocompatible silicone tubing
- Encased radiopaque stripe
- Soft, rounded tip
- Stopcock with latex-free injection site
- 150 ml collection burette with drainage stopcock
- 3 distal openings
- Available in three sizes:  
3.5 French, 5.0 French, 8.0 French

### Ordering Information

#### Urinary Drainage Trays

Model #	Description	Unit of Measure	Image
4193507	3.5 French Uri-Cath™ set with silicone catheter	10/BOX	
4195007	5.0 French Uri-Cath™ set with silicone catheter	10/BOX	
4198007	8.0 French Uri-Cath™ set with silicone catheter	10/BOX	

#### Urinary Drainage Catheters

Model #	Description	Unit of Measure	Image
4193505	3.5 French Uri-Cath™ catheter only	10/BOX	
4195005	5.0 French Uri-Cath™ catheter only	10/BOX	
4198005	8.0 French Uri-Cath™ catheter only	10/BOX	

#### United States:

Utah Medical Products, Inc.  
7043 South 300 West  
Midvale, Utah 84047  
800.533.4984 • 801.566.1200  
www.utahmed.com

#### Europe: EC REP

Utah Medical Products Ltd.  
Athlone Business & Technology Park  
Athlone, Co. Westmeath N37 XK74  
Ireland  
353.90.647.3932



**Gesco™**