



# Nutri-Cath<sup>®</sup> with ENFit<sup>™</sup>

Nutri-Cath with ENFit is a family of locking enteral feeding-only catheters, extension sets, and oral dispensing syringes that:

- \* Complies with ISO 80369-3
- \* Utilizes the well-established GESCO<sup>™</sup> Nutri-Cath specialty enteral feeding catheters
- \* Ensures the avoidance of a possible misconnection with devices administering medications
- \* Ensures a secure connection that minimizes the risk of unwanted disconnections

## 30-Day Catheters

Nutri-Cath's gentle enteral-only catheter is a well-proven solution for either short term or intermittent feedings. Nutri-Cath resists irritation and prevents hyperactive gag reflex during nasogastric, orogastric or nasojejunal placement. Available in silicone and a thermosensitive hi-flow Tecoflex<sup>®</sup> polyurethane, the non-wettable materials will not perforate, clog or adhere to lining of the stomach.

### Silicone Nutri-Cath with ENFit

Packaged 10/bx

|       | 3.5 Fr   | 5.0 Fr   | 6.5 Fr   | 8.0 Fr   |
|-------|----------|----------|----------|----------|
| 40 cm | 4153517E | 4155017E | 4156517E | 4158017E |
| 58 cm | 4153527E | 4155027E | 4156527E | 4158027E |
| 93 cm | 4153537E | 4155037E | 4156537E | 4158037E |

### Polyurethane Nutri-Cath with ENFit

Packaged 10/bx

|       | 3.5 Fr    | 5.0 Fr    | 6.5 Fr    | 8.5Fr     |
|-------|-----------|-----------|-----------|-----------|
| 40 cm | 4153517PE | 4155017PE | 4156517PE | 4158517PE |
| 58 cm | 4153527PE | 4155027PE | 4156527PE | 4158527PE |
| 93 cm | 4153537PE | 4155037PE | 4156537PE | 4158537PE |

## Extension Sets

### Nutri-Cath with ENFit Continuous Feed Extension Sets

Packaged 50/bx

- 4150017E 60" Locking Extension Set, Female ENFit to Male ENFit
- 4150027E 60" Locking Extension Set w/ ENFit side port, Female ENFit to Male ENFit
- 4150067E 30" Locking Extension Set, Female ENFit to Male ENFit
- 415FIT ENFit Transition Connector, Oral Female to Stepped Adaptor



# Oral Dispensing Syringes



Monoject™ dispensing syringes:  
For easy recognition by existing dispensing syringe pumps,  
eliminating potential dosage errors

- \* Purple plungers denote ENFit enteral feeding only
- \* The secure locking ENFit connections are designed to reduce accidental disconnects and misconnects
- \* Eliminates the need for a specialized enteral-only feeding pump

Sterile ENFit oral dispensing syringes, individually packaged

| Volume | Part # |       |
|--------|--------|-------|
| 1 mL:  | 41501E | 60/bx |
| 3 mL:  | 41503E | 80/bx |
| 6 mL:  | 41506E | 80/bx |
| 12 mL: | 41512E | 80/bx |
| 35 mL: | 41535E | 40/bx |
| 60 mL: | 41560E | 30/bx |

Sterile ENFit oral dispensing syringe caps, 5 per package

|     | Part # |         |
|-----|--------|---------|
| CAP | 415CPE | 1000/bx |

ISO 80369-3: 2016 specifies the dimensions and requirements for the design and functional performance of small-bore connectors intended to be used for connections on enteral medical devices and accessories. The ISO 80369-3 connector design in the reverse orientation is commonly known by its federally registered trademarked name ENFit.



UTAH MEDICAL PRODUCTS

FEMCARE



## UNITED STATES

Utah Medical Products, Inc.  
7043 South 300 West  
Midvale, Utah 84047  
Toll-free 800 533 4984  
Tel (801) 566 1200  
Fax (801) 566 2062  
Email [customer@utahmed.com](mailto:customer@utahmed.com)  
[www.utahmed.com](http://www.utahmed.com)

## AUSTRALIA

Femcare Australia  
Unit 12, 5 Gladstone Road  
Castle Hill NSW 2154  
Tel +61 (02) 9045 4100  
Fax +61 (02) 9894 1913  
Email [customer@femcare.com.au](mailto:customer@femcare.com.au)

## EUROPE

Utah Medical Products Ltd.  
Athlone Business & Technology Park  
Garrycastle  
Athlone, Co Westmeath N37 XK74  
Ireland  
Tel +353 (090) 647 3932  
Fax +353 (090) 647 5608  
Email [customerservice\\_ireland@utahmed.com](mailto:customerservice_ireland@utahmed.com)

## Femcare-Nikomed Ltd.

Stuart Court, Spursholt Place, Salisbury Road  
Romsey, Hampshire SO51 6DJ  
United Kingdom  
Tel +44 (01794) 525100  
Fax +44 (01794) 525101  
Email [enquiries@femcare-nikomed.co.uk](mailto:enquiries@femcare-nikomed.co.uk)