

# CORDGUARD®

*The Unified System for Umbilical Cord Clamping,  
Cutting & Blood Collection*

*Designed to  
eliminate  
spraying and  
spattering*

*4-Step Actuation*

*Enclosed collection reservoir  
for obtaining clean,  
uncontaminated blood*

*Ergonomic handle assists  
in delivery of placenta*

*Fully enclosed sharps*

*Reliable, secure  
cord clamp*

*Protective sheath keeps  
tubes clean for labeling*

*Replaceable 5cc vacuum  
tubes conveniently slide  
from handle when full*

*A NEW OBSTETRICAL STANDARD OF CARE*

Umbilical cord management that provides an excellent neonatal blood sample while minimizing exposure to blood.

OSHA  
Regulations  
State:

“All procedures involving blood or other potentially infectious materials shall be performed in such a manner as to minimize splashing, spraying, spattering and generation of droplets of these substances.”

Source: OSHA 29 CFR part 1910.1030  
“Occupational exposure to Bloodborne Pathogens”

United States

Utah Medical Products, Inc.  
7043 South 300 West  
Midvale, Utah 84047  
Toll-free: (800) 533-4984  
Tel: (801) 566-1200  
Fax: (801) 566-2062

Europe

Utah Medical Products Ltd.  
Athlone Business & Technology Park  
Dublin Road  
Athlone, County Westmeath, Ireland  
Tel: (353) 90 647-3932

UTAH MEDICAL  
PRODUCTS INC.

